


Maria Seating System-postural care chair

Two Seating size for selecting

Product size	Small size	Large size
Age range	10 months ~ 7 years old	7 years old ~ 12 years old
Stature range (height)	70cm ~ 120cm	110cm ~ 145cm
Min	70cm	110cm
Max	120cm	145cm
Max user weight	70kgs	70kgs
Seat depth	-1-	
Min	12cm	22cm
Max	35cm	38cm
Max windsweep accommodation	20°	20°
Seat width-with hip laterals	-2-	
Min	22cm	22cm
Max	33cm	37cm
Backrest height	-3-	
Min	26cm	39cm
Max	38cm	45cm
Shoulder width	-4-	
Min	22cm	22cm
Max	38cm	38cm
Chest width	-5-	
Min	20cm	22cm
Max	33cm	33cm
Backrest angle	-6-	90° ~ 150°
Min	13cm	13cm
Max	30cm	36cm
Handle folded up height	-8-	
Min	16cm	16cm
Max	23cm	23cm
Tray size	-9-	31x52cm
Min	16cm	16cm
Max	23cm	23cm
Tray size	-9-	31x52cm
Seat to footplate (without cushion)	-10-	46cm
Min	16cm	16cm
Max	23cm	23cm
Tilt in space	-11-	90° ~ 140°
Min	16cm	16cm
Max	23cm	23cm
Headrest left & right movement	-12-	each of 3cm
Min	16cm	16cm
Max	23cm	23cm
Headrest back & forth movement	-13-	each of 5cm
Min	16cm	16cm
Max	23cm	23cm
Headrest backward angle	-14-	0° ~ 45°

*Seat weight (Seat & Back) 36kg



harvest
LINK INTERNATIONAL

HARVEST EXCEL INTERNATIONAL PTE LTD.
www.harvestlink.com

No. 42, 37th Road, Taichung Industrial Park, Taichung 40768, Taiwan

TEL: +886-4-23580698

FAX: +886-4-23504897

Email: koh6188@ms39.hinet.net / harvest.koh@gmail.com



MARIA Multi-Positioning Chair



MARIA

Multi-Positioning Chair



Upper wide H-shaped chest strap does not affect the smooth of respiratory.



Seat depth can be adjusted from **13~38** cm.
Back height can be adjusted from **39~44** cm.
Foot drop can be adjusted from **13~36** cm.



Caster with lock



Seat/ With tilt-in-space function and can be adjusted different prone.